

## Feature Article

### The Use of Electrolysis for the Treatment of Minor Skin Blemishes

#### From The Editor's Desk

A Ban on Cosmetic Surgery Advertising? Industry says No

#### Latest News

Syneron Buys Ultrashape

Hamilton Fraser Steps in to Rescue Failed Broker Ceart

#### New Product Reviews

Stylage® XXL

#### Food For Thought

Cosmetic Surgery for the Soul? The Role of the Surgeon as Counsellor

#### Legislation

European Standard for Aesthetic Surgery Services - Courting Controversy Already?

#### Eye Spy With My Industry Eye

SMART IDEAS Seminar 2012

#### Business Corner

Facebook Introduces Timelines for Pages

#### Plus our monthly round up

Clinical Information, News Articles & Blogs You May Have Missed and Second Hand Equipment Sales



## New Product Reviews

### Stylage® XXL

STYLAGE® XXL is being called the first 'structural filler' with Hyaluronic Acid (HA).

Aesthetic practitioners now fully understand the important role of facial volume in the ageing process. Volume however is nothing without structure and now Medical Aesthetic Group have launched the first facial restructuring filler or bony structure replacement filler, strong enough yet soft enough to provide structure, lift and longevity.

VIVACY STYLAGE XXL is a highly cross-linked HA with added mannitol designed to provide underlying facial structure in the two most important areas of the face, the cheekbones and chin area, to restore a healthy, balanced face. Mannitol slows down the effects of free radical damage which means the product is not broken down in the hours following the injection and the immediate smoothing and restructuring effect is long-lasting, non-migrating and cohesive.

STYLAGE XXL can be used in combination with all other STYLAGE and HA products. To recreate a full, healthy face they recommend using XXL as a base and then building up volume with XL, L, M or S. The result of the filling can be seen immediately.

STYLAGE XXL can be used for cheek and cheekbone enhancement, chin remodeling, face oval remodeling, mid-face ptosis and facial HIV-associated lipodystrophy treatment.

**For more information please contact: Medical Aesthetic Group, Tel: 02380 676733 [www.magroup.co.uk](http://www.magroup.co.uk)**



### Introducing The New Pixel Q-Switched

Alma Lasers (distributed by ABC Lasers in the UK) have announced that the high powered Pixel Q-Switched is the world's first fractional non ablative 1064nm laser. Used exclusively with the HarmonyXL platform, its unique 2 laser rods provide high power and unmatched treatment efficacy according to the manufacturer.

The high powered Pixel Q-Switched 5x5 mm pixel tip is used for encouraging growth of new collagen fibres, non ablative skin rejuvenation, and is part of the Laser360iQ protocol for complete skin rejuvenation. The Q-Switched laser with non pixel tips is for the treatment of deep pigmentation, tattoo removal, and is part of the Nail360 protocol for the treatment of nail fungus.

Alma Laser's Q-Switched laser module offers tremendous value to practitioners, as it is a single laser module used with a multi application platform. The Q-Switched laser is amendable for use on all parts of the body.

Benefits for the client include: short treatment with minimum down time, unmatched improvement in appearance of wrinkles and skin laxity, fractional non-ablative laser for skin tightening, virtually painless treatment, fewer side effects than CO2 resurfacing and generates new collagen fibres

Benefits for your business include: customisable and upgradeable as your aesthetic practice grows, Plug&Play technology – no downtime when upgrading, no disposables and low cost of ownership – high return on investment and easy to use, safe, versatile and reliable.

**For more information, contact ABC Lasers, email: [info@a-b-c-uk.com](mailto:info@a-b-c-uk.com) web: [www.abclasers.co.uk](http://www.abclasers.co.uk)**

